



# SPRING LOADED T-ADAPTORS

Dynarex Resp-O<sub>2</sub>™ Spring Loaded T-Adaptors are designed to be used in conjunction with a small volume nebulizer (SVN) to help deliver aerosol medication to respiratory vent patients. The T-Adaptor features a spring-loaded valve that opens and closes the nebulizer circuit during insertion and removal.

- Eliminates the need to break the ventilator circuit between treatments
- Ventless design prevents medication leakage
- Accepts most standard aerosol tubing/ventilator circuit
- Not made with natural rubber latex or DEHP

34502

34500

34501

22 mm OD × 22 mm OD

REORDER NO. 34500

PACKAGE QUANTITY  
30/Case

22 mm ID × 22 mm OD

REORDER NO. 34501

PACKAGE QUANTITY  
30/Case

15 mm ID × 15 mm OD

REORDER NO. 34502

PACKAGE QUANTITY  
30/Case